

# Oropharyngeal Diagnostic Swabs


## Description

Oropharyngeal (OP) swabs are used in the oropharynx and anterior nares to collect and carry mucus and saliva in transport tubes to provide a high rate of recovery of microorganisms present, providing accurate analytical test results. US FDA and CDC recommend polyurethane foam swabs for PCR detection of most respiratory viruses including SARS-CoV-2 (COVID-19), rhinoviruses and enteroviruses (Picornaviridae), influenza viruses (Orthomyxoviridae), parainfluenza, metapneumoviruses and respiratory syncytial viruses (Paramyxoviridae), coronaviruses (Coronaviridae) and bacterial infections caused by *Streptococcus pneumoniae* (pneumonia) and *Streptococcus pyogenes* (strep throat).

Foamtec Medical Oropharyngeal swabs utilize Foamtec developed polyurethane foam technology that provides exceptional surface area to capture a maximum quantity of target fluids, while not binding to those fluids providing high subject recovery.

Foamtec Medical proprietary swab cleaning processes and validated sterility ensure no background from RNase and other contaminants for accurate test results. Accurate oropharyngeal swabbing is achieved with rigid, large tipped swabs for accurate placement by the nurse or technician while reaching all areas of the oropharynx and anterior nares.

The plastic handle is engineered to easily snap off so that the swab tip can be deposited into a transport tube.



## Features

- The high surface area captures the maximum amount of mucus and saliva
- No bonding to swabbed body fluids
- Validated sterile, clean swab surfaces
- Rigid handles with large foam tips

## Benefits

- Ability to accurately reach all areas of the oropharynx and anterior nares
- High recovery rates of body fluids for accurate test results
- Maximum patient comfort



**FOAMTEC MEDICAL**  
CONSISTENCY | QUALITY | TRACEABILITY

**Foamtec Medical**  
6575 IH-35 N.,  
Waco, TX 76705


### Foamtec Medical

1-866-922-7160

[www.foamtecmedical.com](http://www.foamtecmedical.com) | [customerservice@foamtecusa.com](mailto:customerservice@foamtecusa.com)


©2023 Foamtec Medical. All rights reserved.

We reserve the right to correct any errors that may occur within this datasheet.




# Oropharyngeal Diagnostic Swabs

## MP1302ST - Packaging Example



## MP1302ST - Swab Tip



## Ordering Information

Part No.	Description	Quantity
MP1302ST	Sterile Oropharyngeal Diagnostic Swab - 1 Swab/Pack, 10 Packs/Box, 100 Boxes/Case	1000 Swabs Total



**FOAMTECMEDICAL**  
CONSISTENCY | QUALITY | TRACEABILITY

**Foamtec Medical**  
6575 IH-35 N.,  
Waco, TX 76705

**Foamtec Medical**  
1-866-922-7160

www.foamtecmedical.com | customerservice@foamtecusa.com