

# Typical Characteristics EZ™ Renew



FOAMTECMEDICAL  
CONSISTENCY | QUALITY | TRACEABILITY

**Embracing sustainability without sacrificing performance. Crafted with renewable materials, offering a sustainable solution in diagnostic applications.**

Type of Foam	Polyether Reticulated Polyurethane
Density	1.80 - 2.20 lb/ft.
Pore Size	75 - 85 ppi
Tear Strength	1.5 lbs/in
Cytotoxicity Grade	0

\*Biocompatible to ISO 10993-5 (non-cytotoxic)

\*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160  
Email: [medicalinquiry@foamtecusa.com](mailto:medicalinquiry@foamtecusa.com)  
Website: [www.foamtecmedical.com](http://www.foamtecmedical.com)

6575 N. IH-35  
Waco, TX 76705  
United States