

Typical Characteristics Vital NP



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

**The softest foam in the Essential
family.**

Type of Foam	Polyether Reticulated Polyurethane
Density	1.25- 1.50 lb/ft.
Pore Size	76 - 95 ppi
Tear Strength	2.0 lbs/in
Cytotoxicity Grade	0

*Biocompatible to ISO 10993-5 (non-cytotoxic)

*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States