

Typical Characteristics EZ™ Clean



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

The perfect balance between comfort and durability, making this foam an ideal choice for daily medical applications. The foam is lightweight and has high performance rates in diagnostic testing.

Type of Foam	Polyester Reticulated Polyurethane
Density	1.80 - 2.20 lb/ft ³
Pore Size	78 - 95 ppi
Tear Strength	3.6 lbs/in
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States