

Typical Characteristics EZ™ Sleek



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

EZ™ Sleek offers a smooth and gentle user experience. This formulation is an ideal choice for diagnostic swabs where patient comfort is valued.

Type of Foam	Polyester Non-Reticulated Polyurethane
Density	4.00 - 5.30 lb/ft ³
Pore Size	75 - 85 ppi
Tear Strength	2.2 lbs/in
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States