

Typical Characteristics EZ™ Test



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

Designed to provide a stable foundation for diagnostic swabs, ensuring accurate and reliable diagnostic medical testing.

Type of Foam	Polyester Reticulated Polyurethane
Density	1.80 - 2.20 lb/ft ³
Pore Size	78 - 95 ppi
Tear Strength	3.6 lbs/in
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States