

Typical Characteristics Echo™



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

Patient comfort can significantly influence the healing process. Echo™ is versatile and can be utilized in positioners for various applications, including imaging, surgery, and rehabilitation. Echo™ foam is commonly used in the orthopedic industry.

Type of Foam	Polyether - Polyurethane
Density	1.30 - 1.50 lb/ft ³
Pore Size	60 - 80 ppi
Tear Strength	N/A
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States