

Typical Characteristics Evo™



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

Gone are the days of cold steel machines for imaging and simple procedures. Thin layers of foam offer a protective layer that increases comfort and minimizes patient
Evo™ is a polyester polyurethane and reticulated foam commonly used in the orthopedic and imaging industry.

Type of Foam	Polyester - Polyurethane
Density	5.40 - 6.60 lb/ft ³
Pore Size	90 ppi
Tear Strength	N/A
Cytotoxicity Grade	2

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States