

Typical Characteristics Imperial™



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

Imperial™ foam is porous and highly dense to provide patient support. It is an excellent option for the personal care environments such as Spa and Wellness Centers, luxury resorts, and private healthcare facilities.

Type of Foam	Polyester Non-Reticulated Polyurethane Foam
Density	3.15 - 3.75 lb/ft ³
Pore Size	85 ppi
Tear Strength	N/A
Cytotoxicity Grade	1

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States