

Typical Characteristics Prima™ X5



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

This foam has the highest density for wound care with a swelling capability.

Type of Foam	Polyether non-Reticulated
Density	4.25-6.25 lb/ft ³
Pore Size	Greater Than 120
Tear Strength	2.5 lbs/in
Cytotoxicity Grade	0

*Biocompatible to ISO 10993-5 (non-cytotoxic)

*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States