

Typical Characteristics UltraPRESS™ S3



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

**Reticulated polyester foam with small pore sizes with great density for
Negative Pressure Wound Therapy.**

Type of Foam	Polyester Reticulated Polyurethane
Density	1.70 - 2.00 lb/ft ³
Pore Size	28 - 38 ppi
Tear Strength	2.5 lbs/in
Cytotoxicity Grade	0

*Biocompatible to ISO 10993-5 (non-cytotoxic)

*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States