

# Typical Characteristics UltraPRESS™ T



FOAMTECMEDICAL  
CONSISTENCY | QUALITY | TRACEABILITY

**Reticulated polyether foam with medium pore size and density that make UltraPRESS™ perfect for Negative Pressure Wound Therapy.**

| Type of Foam       | Polyether Reticulated Polyurethane |
|--------------------|------------------------------------|
| Density            | 1.30 - 1.60 lb/ft <sup>3</sup>     |
| Pore Size          | 40 - 50 ppi                        |
| Tear Strength      | 2.5 lbs/in                         |
| Cytotoxicity Grade | 0                                  |

\*Biocompatible to ISO 10993-5 (non-cytotoxic)

\*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160  
Email: [medicalinquiry@foamtecusa.com](mailto:medicalinquiry@foamtecusa.com)  
Website: [www.foamtecmedical.com](http://www.foamtecmedical.com)

6575 N. IH-35  
Waco, TX 76705  
United States