

Typical Characteristics UltraPRESS™ T3



FOAMTEC MEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

Reticulated polyether foam with medium pore size and a higher density that make UltraPRESS™ perfect for Negative Pressure Wound Therapy

Type of Foam	Polyether Reticulated Polyurethane
Density	1.50 - 1.90 lb/ft ³
Pore Size	20-40 ppi
Cytotoxicity Grade	0

*Biocompatible to ISO 10993-5 (non-cytotoxic)

*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States