

Typical Characteristics UltraSORB® Med



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

We couldn't leave out our UltraSORB line. This foam is non-reticulated polyether and highly absorbent.

Type of Foam	polyether non-reticulated polyurethane
Density	4.00 - 6.00 lb/ft ³
Pore Size	N/A
Tear Strength	3 lbs/in
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States