

Typical Characteristics Vital™ S



FOAMTECMEDICAL
CONSISTENCY | QUALITY | TRACEABILITY

High density and large pores make this foam great for applicators in both everyday life and medical purposes

Type of Foam	Polyether Reticulated Polyurethane
Density	1.80- 2 lb/ft ³
Pore Size	88 - 100 ppi
Tear Strength	2.0 lbs/in
Cytotoxicity Grade	0

***Biocompatible to ISO 10993-5 (non-cytotoxic)**

***Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160
Email: medicalinquiry@foamtecusa.com
Website: www.foamtecmedical.com

6575 N. IH-35
Waco, TX 76705
United States