

# Typical Characteristics Vital™ X



FOAMTECMEDICAL  
CONSISTENCY | QUALITY | TRACEABILITY

**Durable and wicking technology. This foam is easy to cut and shape.**

Type of Foam	Polyether felted Polyurethane
Density	1.25- 1.50 lb/ft <sup>3</sup>
Pore Size	76 - 95 ppi
Tear Strength	2.0 lbs/in
Cytotoxicity Grade	0

**\*Biocompatible to ISO 10993-5 (non-cytotoxic)**

**\*Sterilization - Compatible with Gamma, EO, and steam sterilization (up to 125° C.)**

Phone: +1 866-922-7160  
Email: [medicalinquiry@foamtecusa.com](mailto:medicalinquiry@foamtecusa.com)  
Website: [www.foamtecmedical.com](http://www.foamtecmedical.com)

6575 N. IH-35  
Waco, TX 76705  
United States